

# VERI-Q MULTI



#### • Test parameters



- Auto turn-off
- Auto ejection
- 1000 results memory
- 15  $\mu$ L blood sample (Hb. 7 $\mu$ L)


#### • Accessory



Mini pipette

## Multifunction Monitoring System

- Simple to use for self-test
- Test parameters






CE 0120, MFDS

### Clinical Data




TEST	Lipid			Hemoglobin	Glucose
	Total Cholesterol	Triglyceride	HDL cholesterol		
Measuring range	100~400 mg/dL 2.59~10.4mmol/L	50~600 mg/dL 0.57~6.78mmol/L	20~80 mg/dL 0.52~2.07mmol/L	5~25 g/dL 3.1~15.5mmol/L	20~600 mg/dL 1.11~33.3mmol/L
Measuring time	< 120 Sec.	< 180 Sec.	< 120 Sec.	5 Sec.	5 Sec.
Sample volume	15 µL	15 µL	15 µL	7 µL	0.5 µL
Hematocrit range	30~50%	30~50%	35~50%	-	20~60%
Precision CV (%) Difference (±)	CV <5%	50~100 mg/dL ±10mg/dL 100~600 mg/dL CV <5%	CV <7%	5~10 g/dL ±0.3g/dL 10~25 g/dL CV <4%	<100 mg/dL ±15 mg/dL ≥100 mg/dL CV ≤5%

### Accuracy



Very-Q MULTI shows the acceptance Cobas C111(Roche) for accuracy in clinical evaluation and good performance

### Precision



Very-Q MULTI shows good performance for linearity & repeatability



### Contact Us:

Q-LINE BIOTECH PVT. LTD. (Formerly known as POCT Services Pvt. Ltd.)

Corporate Office:  
20A, Basement, Chhatrapati Shivaji Marg,  
Najafgarh Rd, Near Haldiram,  
New Delhi-110015

+91 11 4557 7407  
sales@qlinebiotech.com  
customercare@qlinebiotech.com  
www.qlinebiotech.com

1800 123 0079

