

# POTROL H 3D

**CONTROL**
**IVD**


Control	Lot	
Low	2602-001L	2026-05-05
Normal	2602-001N	2026-05-05
High	2602-001H	2026-05-05

**Open Vial Stability: 14 days**
**Instrument Name: Medonic M-series/M32**

CONTROL	L
LOT	2602-001L

CONTROL	N
LOT	2602-001N

CONTROL	H
LOT	2602-001H

		Mean	± Range
RBC	10 <sup>12</sup> /L	2.24	0.19
MCV	fL	67.95	5.16
RDW	%	14.66	6.61
Hct	%	15.20	1.70
Plt	10 <sup>9</sup> /L	60.79	25.20
MPV	fL	8.58	1.89
WBC	10 <sup>9</sup> /L	1.82	0.60
Hgb	g/dL	5.96	0.56
MCH	pg	26.73	1.96
MCHC	g/dL	39.37	3.87
Lymph#	10 <sup>9</sup> /L	1.06	0.50
Mid#	10 <sup>9</sup> /L	0.16	0.25
Gran#	10 <sup>9</sup> /L	0.60	0.20
Lymph	%	55.44	12.96
Mid	%	8.92	5.18
Gran	%	35.64	13.60

	Mean	± Range
	4.56	0.31
	76.81	5.64
	14.33	7.44
	35.07	4.03
	223.53	52.76
	8.56	1.96
	7.34	1.28
	13.15	1.17
	28.85	1.95
	37.52	3.57
	1.63	0.82
	0.75	0.41
	4.95	0.92
	22.66	8.08
	10.46	4.86
	66.89	8.89

	Mean	± Range
	5.63	0.48
	85.33	7.72
	14.54	7.66
	48.04	7.05
	532.82	88.06
	8.94	4.25
	20.99	3.15
	17.96	2.02
	31.91	1.90
	37.42	3.87
	2.61	0.84
	1.07	0.58
	17.28	2.72
	12.65	3.23
	5.24	2.44
	82.11	4.73


**Q-Line Biotech Limited.**

 Plot No. 3105, Gindan Khera,  
 Amausi Station Road, Nadarganj,  
 Lucknow, Uttar Pradesh-226008,  
 India.

Prepared By:

  
 13/02/26

Approved By:

  
 13/02/26

*Note: The assay sheet is available for download from the official Q-Line website at <https://qlinebiotech.com/documents>.*